

# Anti-S100A12 antibody



Catalog Number: 100444

## Product name

Anti-S100A12 antibody

## Immunogen

[Human S100A12 recombinant protein](#)

## Specificity

Human S100A12

## Antibody description

Rabbit monoclonal to S100A12

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human S100A12 (rh S100A12; NP\_005612.1; Met1-Glu92).

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 2-10  $\mu$ g/mL

## Validations



S100A12 Antibody, Rabbit MAb,  
Immunohistochemistry

Immunochemical staining of human S100A12 in human tonsil with rabbit monoclonal antibody (5  $\mu$ g/mL, formalin-fixed paraffin embedded sections).