

Anti-Layilin / LAYN antibody



Catalog Number: 100567

Product name

Anti-Layilin / LAYN antibody

Immunogen

[Rat Layilin / LAYN \(His Tag\) recombinant protein](#)

Specificity

Rat Layilin / LAYN

Antibody description

Rabbit polyclonal to Layilin / LAYN

Preparation

Produced in rabbits immunized with purified, recombinant Rat Layilin / LAYN (rh Layilin / LAYN; NP_001178926.1; Met1-Glu224). Layilin / LAYN specific IgG was purified by Rat Layilin / LAYN affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA:0.1-0.2 μ g/mL

This antibody can be used at 0.1-0.2 μ g/mL with the appropriate secondary reagents to detect Rat Layilin / LAYN. The detection limit for Rat Layilin / LAYN is < 0.039 ng/well.

IHC-P: 0.1-2 μ g/ml

Validations



Layilin / LAYN Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry