

# Anti-Layilin / LAYN antibody



Catalog Number: 100567

## Product name

Anti-Layilin / LAYN antibody

## Immunogen

[Rat Layilin / LAYN \(His Tag\) recombinant protein](#)

## Specificity

Rat Layilin / LAYN

## Antibody description

Rabbit polyclonal to Layilin / LAYN

## Preparation

Produced in rabbits immunized with purified, recombinant Rat Layilin / LAYN (rh Layilin / LAYN; NP\_001178926.1; Met1-Glu224). Layilin / LAYN specific IgG was purified by Rat Layilin / LAYN affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC-P


## Dilutions

ELISA:0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat Layilin / LAYN. The detection limit for Rat Layilin / LAYN is < 0.039 ng/well.

IHC-P: 0.1-2 µg/ml

## Validations



Layilin / LAYN Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry