

# Anti-GATA3 antibody



Catalog Number: 100569

## Product name

Anti-GATA3 antibody

## Specificity

Human GATA3

## Antibody description

Rabbit polyclonal to GATA3

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human GATA3, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB10-20 µg/ml

## Validations



Lanes	A	B	C
Items			
Sample (whole cell lysate)	Jurkat	HeLa	MCF7
Sample Volume (µg/lane)	30	30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	10-20 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+), at 1:5000 dilution.		
Developed using Odyssey imaging system.			
Explanation	Predicted band size : 48 kDa Observed band size : 45 kDa		

GATA3 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot