

# Anti-Annexin A8 / ANXA8 antibody



Catalog Number: 100735

## Product name

Anti-Annexin A8 / ANXA8 antibody

## Immunogen

[Human Annexin A8 / ANXA8 \(His Tag\) recombinant protein](#)

## Specificity

Human Annexin A8 / ANXA8

## Antibody description

Rabbit polyclonal to Annexin A8 / ANXA8

## Preparation

Produced in rabbits immunized with purified, recombinant Human Annexin A8 / ANXA8 (rh Annexin A8 / ANXA8; AAH73755.1; Met1-Pro327). Annexin A8 / ANXA8 specific IgG was purified by Human Annexin A8 / ANXA8 affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB


## Dilutions

WB: 10-20 µg/ml

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human Annexin A8 / ANXA8. The detection limit for Human Annexin A8 / ANXA8 is < 0.039 ng/well.

## Validations



Lanes	A
Items	
Sample (whole cell lysate)	Junkat
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	10-20 µg/ml
Secondary Antibody	Dylight 800 labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution
Developed using	Odyssey imaging system
Explanation	Predicted band size : 37 kDa Observed band size : 37 kDa

Annexin A8 / ANXA8 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot