

Anti-BCAM / CD239 antibody



Catalog Number: 100784

Product name

Polyclonal

Anti-BCAM / CD239 antibody

Ig Type

Immunogen

Rabbit IgG

[Rat BCAM / CD239 \(His Tag\) recombinant protein](#)

Applications

ELISA, IHC-P

Specificity

Rat BCAM

Dilutions

ELISA: 0.1-0.2 µg/mL

Antibody description

Rabbit polyclonal to BCAM / CD239

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat BCAM. The detection limit for Rat BCAM is approximately 0.00245 ng/well.

Preparation

Produced in rabbits immunized with purified, recombinant Rat BCAM (rR BCAM; Q9ESS6; Met1-Ala543). BCAM specific IgG was purified by Rat BCAM affinity chromatography.

IHC-P: 0.1-2 µg/mL

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Validations


Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality



BCAM Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of rat BCAM in rat kidney with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).