

# Anti-FABP3/H-FABP antibody



Catalog Number: 100794

## Product name

Anti-FABP3/H-FABP antibody

## Specificity

Human FABP3 / H-FABP

## Antibody description

Rabbit polyclonal to FABP3/H-FABP

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human FABP3 / H-FABP, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 0.5-5 µg/ml


IP: 0.5-2 uL/mg of lysate

## Validations



Items	Lanes
Sample	A Mouse heart tissue lysate
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	0.5-5 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution. Developed using Odyssey imaging system.
Explanation	Predicted band size : 15 kDa Observed band size : 13 kDa

## FABP3 / H-FABP Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



Items	Lanes
Sample (whole cell lysate)	A Jurkat
Sample quantity	0.5 mg
IP antibody quantity	2 µg
Immunomagnetic beads (Protein G)	60 mg
Gel	13% SDS-PAGE reducing gel
Primary antibody	KL1F-SMCC-62N0A751 antibody at 10 µg/ml
Secondary antibody	Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:2000 dilution

## FABP3 / H-FABP Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation