

Anti-ST6GAL1/ST6GAL-I antibody



Catalog Number: 100800

Product name

Anti-ST6GAL1/ST6GAL-I antibody

Immunogen

[Human ST6GAL1/ST6GAL-I \(His Tag\) recombinant protein](#)

Specificity

Human ST6GAL1 / ST6GAL-I

Antibody description

Rabbit monoclonal to ST6GAL1/ST6GAL-I

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human ST6GAL1 / ST6GAL-I (rh ST6GAL1 / ST6GAL-I; P15907; Lys27-Cys406).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10 µg/mL

Validations



ST6GAL1 / CD75 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human ST6GAL1 in human kidney with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).

ST6GAL1 / CD75 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human ST6GAL1 in human prostate with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to epithelium.