

Anti-CD98/SLC3A2 antibody



Catalog Number: 100882

Product name

Anti-CD98/SLC3A2 antibody

Immunogen

[Mouse CD98/SLC3A2 \(His Tag\) recombinant protein](#)

Specificity

Mouse SLC3A2

No cross-reactivity in ELISA with Human SLC3A2

Antibody description

Rabbit monoclonal to CD98/SLC3A2

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse SLC3A2 (rM SLC3A2; NP_001154885.1; Ala139-Ala565).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 1-10 µg/mL

Validations



SLC3A2 / CD98 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse SLC3A2 in mouse kidney with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to nephric tubule.