

# Anti-Podoplanin antibody



Catalog Number: 100930

## Product name

Anti-Podoplanin antibody

## Immunogen

[Mouse Podoplanin \(His Tag\) recombinant protein](#)

## Specificity

Mouse Podoplanin / PDPN

## Antibody description

Rabbit monoclonal to Podoplanin

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Podoplanin / PDPN (rM Podoplanin / PDPN; NP\_034459.2; Met1-Lys133).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 1-10 µg/mL

## Validations



Podoplanin / PDPN Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse PDPN in mouse brain with rabbit monoclonal antibody (2 µg/mL, formalin-fixed paraffin embedded



sections).

Podoplanin / PDPN Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse PDPN in mouse kidney with rabbit monoclonal antibody (6 µg/mL, formalin-fixed paraffin embedded sections).