

# Anti-JAG1/Jagged 1(CD339) antibody



Catalog Number: 101051

## Product name

Anti-JAG1/Jagged 1(CD339) antibody

## Immunogen

[Human JAG1/Jagged 1\(CD339\) \(His Tag\) recombinant protein](#)

## Specificity

Human JAG1 / Jagged 1

## Antibody description

Mouse monoclonal to JAG1/Jagged 1(CD339)

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human JAG1 / Jagged 1 (rh JAG1 / Jagged 1; NP\_000205.1; Met1-Ser1046) and conjugated with PE under optimum conditions, the unreacted PE was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

FCM

## Dilutions

FCM

## Validations



JAG1 / Jagged 1 Antibody (PE), Mouse MAb, Flow cytometric

Profile of anti-JAG1(CD339) reactivity on HepG2 cells analyzed by flow cytometry.