

# Anti-Thrombomodulin / THBD antibody



Catalog Number: 101247

## Product name

Rabbit IgG

Anti-Thrombomodulin / THBD antibody

## Applications

IHC-P

## Immunogen

[Rat Thrombomodulin / THBD \(His Tag\) recombinant protein](#)

## Dilutions

IHC-P: 0.5-5 µg/mL


## Specificity

Rat Thrombomodulin / THBD

## Validations

## Antibody description

Rabbit monoclonal to Thrombomodulin / THBD



## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat Thrombomodulin / THBD (rR Thrombomodulin / THBD; O35370; Met1-Ser517).

Thrombomodulin / THBD Antibody, Rabbit MAb, Immunohistochemistry

## Formulation

0.2 µm filtered solution in PBS


Immunochemical staining of rat THBD in rat lung with rabbit monoclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.



Thrombomodulin / THBD Antibody, Rabbit MAb, Immunohistochemistry

## Clonality

Monoclonal

Immunochemical staining of rat THBD in rat liver with rabbit monoclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).

## Ig Type