

Anti-CD83 antibody



Catalog Number: 101254

Product name

Anti-CD83 antibody

Immunogen

[Human CD83 \(His Tag\) recombinant protein](#)

Specificity

Human CD83

Antibody description

Rabbit monoclonal to CD83

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CD83 (rh CD83; Q01151-1; Met1-Ala143).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-10 µg/ml

Validations



Lanes	A	B
Items		
Sample (whole cell lysate)	Jurkat	293T
Sample Volume (µg/lane)	30	30
Gel	13% SDS-PAGE reducing gel	
Recommended Concentration	5-10 µg/ml	
Secondary Antibody	DyLight 800 labeled Antibody to Rabbit IgG (H+L) at 1:5000 dilution.	
Developed using Odyssey imaging system.		
Explanation	Predicted band size : 23 kDa Observed band size : 55 kDa	

CD83 Antibody, Rabbit MAb, Western blot