

## Product name

Anti-TNFR1/CD120a/TNFRSF1A antibody

## Immunogen

[Human TNFR1/CD120a/TNFRSF1A \(His Tag\) recombinant protein](#)

## Specificity

Human TNFR1 / CD120a / TNFRSF1A  
**No cross-reactivity** with Mouse TNFR1 (504575) in ELISA assay

## Antibody description

Rabbit monoclonal to TNFR1/CD120a/TNFRSF1A

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human TNFR1 / CD120a / TNFRSF1A (rh TNFR1 / CD120a / TNFRSF1A; NP\_001056.1; Met1-Thr211) and was produced using recombinant antibody technology.

## Formulation

0.2 µm filtered solution in Histidine and Arginine buffer containing 120mM NaCl, 0.02% Tween 80, pH6.0

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


B/N

## Dilutions

**Block** - This TNFR1 neutralizing antibody can block the binding of recombinant TNFa to its recombinant receptor TNFR1.

**Neutralization** - The neutralization activity of TNFR1/TNFRSF1A Neutralizing Antibody is Measured by its ability to neutralize TNFR1/TNFRSF1A-mediated inhibition of cytotoxicity in the L-929 mouse fibroblast cell line. The Neutralization titer (IC50) is typically 0.5-1.5 µg/mL in the presence of 0.3 µg/mL Recombinant Human TNFR1/TNFRSF1A, 0.2 ng/mL Recombinant Human TNFa and 1 µg/mL actinomycin D.

## Validations



TNFR1 / CD120a / TNFRSF1A Neutralizing Antibody

**TNFR1/TNFRSF1A-mediated inhibition of cytotoxicity was Neutralized by Human TNFR1 Antibody.** Recombinant Human TNFR1/TNFRSF1A (Catalog 10872-H08H) inhibits Recombinant Human TNFa (Catalog 10602-HNAE) induced cytotoxicity in the L-929 mouse fibroblast cell line. Inhibition of Recombinant Human TNFa (0.2 ng/mL) activity elicited by

# Anti-TNFR1/CD120a/TNFRSF1A antibody



Catalog Number: 101281

---

Recombinant Human TNFR1/TNFRSF1A (0.3 µg/mL) is neutralized by increasing concentrations of Human TNFR1/TNFRSF1A Monoclonal

Antibody (Catalog 10872-R111). The IC<sub>50</sub> is typically 0.5-1.5 µg/mL in the presence of the metabolic inhibitor actinomycin D (1 µg/mL).