

# Anti-B2M/beta-2 microglobulin antibody



Catalog Number: 101304

## Product name

Anti-B2M/beta-2 microglobulin antibody

## Immunogen

[Cynomolgus B2M/beta-2 microglobulin \(His Tag\) recombinant protein](#)

## Specificity

Cynomolgus B2M / beta-2 microglobulin

## Antibody description

Rabbit polyclonal to B2M/beta-2 microglobulin

## Preparation

Produced in rabbits immunized with purified, recombinant Cynomolgus B2M / beta-2 microglobulin (Q6V7J5; Met1-Met119). B2M / beta-2 microglobulin specific IgG was purified by Cynomolgus B2M / beta-2 microglobulin affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions


WB: 2-10 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Cynomolgus B2M / beta-2 microglobulin. The detection limit for Cynomolgus B2M / beta-2 microglobulin is approximately < 0.039 ng/well.


IP: 1-4 µg/mg of lysate

## Validations



Items	Lanes	A	B	C	D
Sample (whole cell lysate)		A431	NIH-3T3	HeLa	HL60
Sample Volume (µg/lane)		30	30	30	30
Gel	13% SDS-PAGE reducing gel				
Recommended Concentration	2-10 µg/ml				
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.				
Explanation	Developed using Odyssey imaging system. Predicted band size : 14 kDa Observed band size : 14 kDa				

B2M / beta-2 microglobulin Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



Items	Lanes	A
Sample (whole cell lysate)		Jurkat
Sample quantity		0.5 mg
IP antibody quantity		2 µg
Protein G agarose		15 µl of 50% Protein G Agarose
Gel		13% SDS-PAGE reducing gel
Primary antibody		Cynomolgus B2M antibody at 10 µg/ml
Secondary antibody		Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution.

B2M / beta-2 microglobulin Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation

Immunochemical staining of human CLSTN1 in

# Anti-B2M/beta-2 microglobulin antibody



Catalog Number: 101304

---

human brain with rabbit polyclonal antibody (0.2  $\mu\text{g}/\text{mL}$ , formalin-fixed paraffin embedded sections).