

Anti-PERP1 / MZB1 antibody



Catalog Number: 101671

Product name

Anti-PERP1 / MZB1 antibody

Immunogen

[Human PERP1 / MZB1 \(His Tag\) recombinant protein](#)

Specificity

Human PERP1 / MZB1

Antibody description

Rabbit monoclonal to PERP1 / MZB1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human PERP1 / MZB1 (rh PERP1 / MZB1; Q8WU39-1; Met 1- Thr 185).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10 µg/mL

Validations



PERP1 / MZB1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human MGC29506 in human spleen with rabbit monoclonal antibody (2 µg/mL, formalin-fixed paraffin embedded



sections).

PERP1 / MZB1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human MGC29506 in human small intestine with rabbit monoclonal antibody (2 µg/mL, formalin-fixed paraffin embedded sections).