

# Anti-Neuroigin 1/NLGN1 antibody



Catalog Number: 101707

## Product name

Anti-Neuroigin 1/NLGN1 antibody

## Immunogen

[Human Neuroigin 1/NLGN1 \(His Tag\) recombinant protein](#)

## Specificity

Human Neuroigin 1 / NLGN1

## Antibody description

Rabbit monoclonal to Neuroigin 1/NLGN1

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Neuroigin 1 / NLGN1 (rh Neuroigin 1 / NLGN1; NP\_055747.1; Met1-Ser677).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-10 µg/mL

## Validations



Neuroigin 1 / NLGN1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human NLGN1 in human brain with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).