

# Anti-DMBT1/Muclin antibody



Catalog Number: 102735

## Product name

Anti-DMBT1/Muclin antibody

## Immunogen

[Human DMBT1/Muclin \(His Tag\) recombinant protein](#)

## Specificity

Human DMBT1 / GP340

## Antibody description

Rabbit polyclonal to DMBT1/Muclin

## Preparation

Produced in rabbits immunized with purified, recombinant Human DMBT1 (rh DMBT1; NP\_004397.2; Met 1-Ser 220). DMBT1 specific IgG was purified by human DMBT1 affinity chromatography

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC-P


## Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human DMBT1. The detection limit for Human DMBT1 is approximately 0.00975 ng/well.

IHC-P: 0.1-2 µg/mL

## Validations



DMBT1 / Muclin Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human DMBT1 in human colon carcinoma with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).