

Anti-GFRA1 antibody



Catalog Number: 103192

Product name

Anti-GFRA1 antibody

Immunogen

[Mouse GFRA1 \(His Tag\) recombinant protein](#)

Specificity

Mouse GFRA1

Antibody description

Rabbit monoclonal to GFRA1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse GFRA1 (rM GFRA1; AAB86600.1; Met 1-Ser 425).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse GFRA1. The detection limit for Mouse GFRA1 is approximately 0.0049 ng/well.


IHC-P: 1-10 µg/mL

Validations



GFRA1 Antibody, Rabbit MAb, Immunochemistry

Immunochemical staining of mouse GFRA1 in mouse brain with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).

GFRA1 Antibody, Rabbit MAb, Immunochemistry

Immunochemical staining of mouse GFRA1 in mouse kidney with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).