

Anti-CXCL9 / MIG antibody



Catalog Number: 103714

Product name

Anti-CXCL9 / MIG antibody

Immunogen

[Human CXCL9 / MIG recombinant protein](#)

Specificity

Human CXCL9 / MIG

Antibody description

Rabbit polyclonal to CXCL9 / MIG

Preparation

Produced in rabbits immunized with purified, recombinant Human CXCL9 / MIG (rh CXCL9 / MIG; Q07325; Thr 23-Thr 125). CXCL9 / MIG specific IgG was purified by Human CXCL9 / MIG affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 5-10 μ g/mL

ELISA: 0.1-0.2 μ g/mL

This antibody can be used at 0.1-0.2 μ g/mL with the appropriate secondary reagents to detect Human CXCL9. The detection limit for Human CXCL9 is approximately 0.00975 ng/well.

Validations



Items	Lanes	A
Sample (whole cell lysate)		A431
Sample Volume (μ g/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		5-10 μ g/ml
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.

CXCL9 / MIG Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot