

Anti-S100A8 antibody



Catalog Number: 103809

Product name

Anti-S100A8 antibody

Clonality

Monoclonal

Immunogen

[Human S100A8 \(His Tag\) recombinant protein](#)

Ig Type

Rabbit IgG

Specificity

Human S100A8

Applications

IHC-P

Antibody description

Rabbit monoclonal to S100A8


Dilutions

IHC-P: 1-10 µg/mL

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human S100A8 (rh S100A8; NP_002955.2; Met 1-Glu 93).

Validations



Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

S100A8 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of human S100A8 in human spleen with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).