

Anti-IL-1R1/CD121a antibody



Catalog Number: 104054

Product name

Anti-IL-1R1/CD121a antibody

Immunogen

[Rat IL-1R1/CD121a \(His Tag\) recombinant protein](#)

Specificity

Rat IL-1R1 / CD121a

Antibody description

Rabbit monoclonal to IL-1R1/CD121a

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat IL-1R1 / CD121a (rR IL-1R1 / CD121a; NP_037255.3; Met1-Lys352) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



IL-1R1 / CD121a Antibody (FITC), Rabbit MAb, Flow cytometric analysis

Flow cytometric analysis of Rat IL1R1(CD121a) expression on rat splenocytes. SD splenocytes were stained with FITC-conjugated anti-Rat IL1R1(CD121a). The histogram were derived from gated events with the forward and side light-scatter characteristics of intact cells.