

Anti-IL-1R1/CD121a antibody



Catalog Number: 104168

Product name

Anti-IL-1R1/CD121a antibody

Immunogen

[Mouse IL-1R1/CD121a \(His Tag\) recombinant protein](#)

Specificity

Mouse IL-1R1 / CD121a

Antibody description

Rabbit monoclonal to IL-1R1/CD121a

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL-1R1 / CD121a (rM IL-1R1 / CD121a; P13504; Met1-Lys338).

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM: 0.5-2 μ g/Test

Validations



IL-1R1 / CD121a Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of IL1R1(CD121a) expression on NIH3T3 cells. Cells were stained with purified anti-Mouse IL1R1(CD121a), then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.