

# Anti-RBP4 antibody



Catalog Number: 104607

## Product name

Anti-RBP4 antibody

## Immunogen

[Mouse RBP4 \(His Tag\) recombinant protein](#)

## Specificity

Mouse RBP4

## Antibody description

Rabbit monoclonal to RBP4

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse RBP4 (rM RBP4; NP\_035385.1; Met1-Leu201) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM

## Validations



RBP4 Antibody (FITC), Rabbit MAb, Flow cytometric

Flow cytometric analysis of Mouse RBP4 expression in spleen lymphocytes. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and then stained with FITC Rabbit anti-Mouse RBP4. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.