

Anti-Contactin 3/CNTN3 antibody



Catalog Number: 105321

Product name

Anti-Contactin 3/CNTN3 antibody

Immunogen

[Rat Contactin 3/CNTN3 \(His Tag\) recombinant protein](#)

Specificity

Rat Contactin 3 / CNTN3

Antibody description

Rabbit polyclonal to Contactin 3/CNTN3

Preparation

Produced in rabbits immunized with purified, recombinant Rat Contactin 3 / CNTN3 (rh Contactin 3 / CNTN3; Q62682; Met1-Gly1001). Contactin 3 / CNTN3 specific IgG was purified by Rat Contactin 3 / CNTN3 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat Contactin 3 / CNTN3. The detection limit for Rat Contactin 3 / CNTN3 is < 0.039 ng/well.

IHC-P: 0.1-2 µg/mL

Validations



Contactin 3 / CNTN3 Antibody, Rabbit PAb, Antigen Affinity Purified

Immunochemical staining of rat CNTN3 in rat brain with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).