

General Information

Product name:

Anti-CaM Kinase 4/CaMKIV antibody

Immunogen:

[Mouse CaM Kinase 4/CaMKIV \(His & GST Tag\) recombinant protein](#)

Applications:

ELISA

Specificity:

Mouse CaM Kinase 4 / CaMKIV

Antibody description:

Rabbit monoclonal to CaM Kinase 4/CaMKIV

Properties

Preparation:

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CaM Kinase 4 / CaMKIV (rh CaM Kinase 4 / CaMKIV; P08414; Met1-Tyr469).

Formulation:

0.2 µm filtered solution in PBS

Storage:

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality:

Monoclonal

Ig Type:Rabbit IgG

Applications

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse CaM Kinase 4 / CaMKIV. The detection limit for Mouse CaM Kinase 4 / CaMKIV is 0.0195 ng/well.

Validations