

Anti-SIGIRR/TIR8 antibody



Catalog Number: 105848

Product name

Anti-SIGIRR/TIR8 antibody

Immunogen

[Human SIGIRR/TIR8 \(His Tag\) recombinant protein](#)

Specificity

Human SIGIRR / TIR8

Antibody description

Rabbit monoclonal to SIGIRR/TIR8

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human SIGIRR / TIR8 (rh SIGIRR / TIR8; Q6IA17-1; Met1-His118) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



SIGIRR / TIR8 Antibody (PE), Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Human SIGIRR expression on SW480 cells. Cells were stained with PE-conjugated anti-Human SIGIRR. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.