

Anti-EPHB3 antibody



Catalog Number: 105877

Product name

Anti-EPHB3 antibody

Immunogen

[Mouse EphB3/HEK2/Eph Receptor B3 \(His Tag\) recombinant protein](#)

Specificity

Mouse EphB3 / HEK2

Antibody description

Rabbit Monoclonal to Mouse EPHB3

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse EphB3 / HEK2 (rM EphB3 / HEK2; NP_034273.1; Met1-Thr537) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



EphB3 / HEK2 Antibody (PE), Rabbit MAb, Flow cytometric

Flow cytometric analysis of Mouse EphB3 expression in CoLo205 cells.