

Anti-IL20RB antibody



Catalog Number: 106024

Product name

Anti-IL20RB antibody

Immunogen

[Rat IL20RB \(His Tag\) recombinant protein](#)

Specificity

Rat IL20RB

Antibody description

Rabbit polyclonal to IL20RB

Preparation

Produced in rabbits immunized with purified, recombinant Rat IL20RB (rh IL20RB; B8K1X9; Met1-Pro227). IL20RB specific IgG was purified by Rat IL20RB affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 10-20 µg/ml

ELISA:0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat IL20RB. The detection limit for Rat IL20RB is 0.00975 ng/well.

Validations



Items	Lanes	A
Sample (whole cell lysate)		HeLa
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10-20 µg/ml
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.

IL20RB Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot