

Anti-EphB3/HEK2/Eph Receptor B3 antibody



Catalog Number: 106061

Product name

Anti-EphB3/HEK2/Eph Receptor B3 antibody

Immunogen

[Mouse EphB3/HEK2/Eph Receptor B3 \(His Tag\) recombinant protein](#)

Specificity

Mouse EphB3 / HEK2

Antibody description

Rabbit monoclonal to EphB3/HEK2/Eph Receptor B3

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse EphB3 / HEK2 (rM EphB3 / HEK2; NP_034273.1; Met1-Thr537) and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of

water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



EphB3 / HEK2 Antibody (APC), Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Mouse EphB3 expression on CoLo205 cells. Cells were stained with APC-conjugated anti-Mouse EphB3. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.