

Anti-CD112/Nectin-2/PVRL2 antibody



Catalog Number: 106118

Product name

Anti-CD112/Nectin-2/PVRL2 antibody

Immunogen

[Rat CD112/Nectin-2/PVRL2 \(His Tag\) recombinant protein](#)

Specificity

Rat CD112 / Nectin-2 / PVRL2

Antibody description

Rabbit polyclonal to CD112/Nectin-2/PVRL2

Preparation

Produced in rabbits immunized with purified, recombinant Rat CD112 / Nectin-2 / PVRL2 (rh CD112 / Nectin-2 / PVRL2; Q9ESS6; Met1-Gly351). CD112 / Nectin-2 / PVRL2 specific IgG was purified by Rat CD112 / Nectin-2 / PVRL2 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA:0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat CD112 / Nectin-2 / PVRL2. The detection limit for Rat CD112 / Nectin-2 / PVRL2 is < 0.039 ng/well.

IHC-P: 1-5 µg/ml

Validations



CD112 / Nectin-2 / PVRL2 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry