

# Anti-SFTPD / SP-D / SFTP4 antibody



Catalog Number: 106195

## Product name

Anti-SFTPD / SP-D / SFTP4 antibody

## Immunogen

[Mouse SFTPD / SP-D / SFTP4 \(His Tag\) recombinant protein](#)

## Specificity

Mouse SFTPD / SP-D / SFTP4

## Antibody description

Rabbit monoclonal to SFTPD / SP-D / SFTP4

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse SFTPD / SP-D / SFTP4 (rM SFTPD / SP-D / SFTP4; NP\_033186.1; Met1-Phe374).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.5-5 µg/mL

## Validations



SFTPD / SP-D / SFTP4 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse SFTPD in mouse lung with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).