

# Anti-PIGR antibody



Catalog Number: 106209

## Product name

Anti-PIGR antibody

## Immunogen

[Human PIGR \(365 Ser/Gly, His Tag\) recombinant protein](#)

## Specificity

Human PIGR

## Antibody description

Rabbit monoclonal to PIGR

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human PIGR (rh PIGR; NP\_002635.2; Met1-Arg638,365Ser/Gly).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 15-30 µg/mL

## Validations



PIGR Antibody, Rabbit MAb, Immunochemistry

Immunochemical staining of human PIGR in human liver (cirrhosis) with rabbit monoclonal antibody (15 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to hepatocyte.



PIGR Antibody, Rabbit MAb, Immunochemistry

Immunochemical staining of human PIGR in human breast carcinoma with rabbit monoclonal antibody (15 µg/mL, formalin-fixed paraffin embedded sections).