

Anti-ASGPR1 antibody



Catalog Number: 106339

Product name

Anti-ASGPR1 antibody

Immunogen

[Mouse ASGPR1 \(His Tag\) recombinant protein](#)

Specificity

Mouse ASGPR1 / ASGR1

Antibody description

Rabbit monoclonal to ASGPR1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse ASGPR1 / ASGR1 (rM ASGPR1 / ASGR1; NP_033844.1; Ser60-Asn284).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.5-2 µg/mL

Validations



ASGPR1 / ASGR1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse ASGR1 in mouse liver with rabbit monoclonal antibody (2 µg/mL, formalin-fixed paraffin embedded sections).