

Anti-UCHL3 antibody



Catalog Number: 106400

Product name

Anti-UCHL3 antibody

Immunogen

[Mouse UCHL3 \(His Tag\) recombinant protein](#)

Specificity

Mouse UCHL3

No cross-reactivity in ELISA with Human UCHL3

Antibody description

Rabbit monoclonal to UCHL3

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse UCHL3 (rM UCHL3; Q9JKB1; Glu2-Ala230).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10 µg/mL

Validations



UCHL3 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of mouse UCHL3 in mouse bone with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).


UCHL3 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of mouse UCHL3 in mouse brain with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded

Anti-UCHL3 antibody



Catalog Number: 106400



UCHL3 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of mouse UCHL3 in
mouse heart with rabbit monoclonal antibody (5
 $\mu\text{g}/\text{mL}$, formalin-fixed paraffin embedded
sections).

sections).