

Anti-CD38 antibody



Catalog Number: 106426

General Information

Product name:

Anti-CD38 antibody

Immunogen:

[Human CD38 \(His Tag\) recombinant protein](#)

Applications:

IHC-P

Specificity:

Human CD38

No cross-reactivity in ELISA with Mouse CD38

Antibody description:

Rabbit monoclonal to CD38

Properties

Preparation:

This antibody was obtained from a rabbit immunized with purified, recombinant Human CD38 (rh CD38; NP_001766.2; Val43-Ile300).

Formulation:

0.2 µm filtered solution in PBS with 5% trehalose

Storage:

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody

products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality:


Monoclonal

Ig Type:Rabbit IgG

Applications

IHC-P: 1-10 µg/mL

Validations



CD38 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of human CD38 in human tonsil with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).