

# Anti-Transferrin/TF antibody



Catalog Number: 106510

## Product name

Anti-Transferrin/TF antibody

## Immunogen

[Mouse Transferrin/TF \(His Tag\) recombinant protein](#)

## Specificity

Mouse Transferrin / TF

## Antibody description

Rabbit polyclonal to Transferrin/TF

## Preparation

Produced in rabbits immunized with purified, recombinant Mouse Transferrin / TF (rh Transferrin / TF; Q921I1; Met1-His697). Transferrin / TF specific IgG was purified by Mouse Transferrin / TF affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.1-2 µg/ml

## Validations



Transferrin / TF Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry