

# CAPZA1 antibody



Catalog Number: 107032

## Product name

CAPZA1 antibody

## Immunogen

[Human CAPZA1 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

CAPZA1 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Positive IF detected in MCF-7 cells, Hela cells. Observed molecular weight by Western-blot: 37kD

## Preparation

This antibody was obtained by immunization of CAPZA1 recombinant protein (Accession Number: NM\_006135). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107032(CAPZA1 antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells, using CAPZA1 antibody Catalog No: at 1:50 dilution and Rhodamine-labeled goat anti-mouse IgG (red).