

CD22 antibody



Catalog Number: 107053

Product name

CD22 antibody

WB: 1:500-1:5000

IHC: 1:100-1:800

Immunogen

Specificity

Human; other species not tested.

Antibody description

CD22 Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 135 kDa

Preparation

This antibody was obtained by immunization of CD22 recombinant protein (Accession Number: NM_001278417). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107053(CD22 antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107053(CD22 Antibody) at dilution of 1:400 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



CD22 antibody



Catalog Number: 107053

Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107053(CD22 Antibody) at dilution of 1:400

(under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).