

CD22 antibody



Catalog Number: 107053

Product name

CD22 antibody

IHC: 1:100-1:800


Specificity

Human; other species not tested.

Validations

Antibody description

CD22 Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 135 kDa




Preparation

This antibody was obtained by immunization of CD22 recombinant protein (Accession Number: NM_001278417). Purification method: Protein A purified.

Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107053(CD22 antibody) at dilution of 1:1500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Monoclonal

Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107053(CD22 Antibody) at dilution of 1:400 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).

Ig Type

Mouse IgG2a



Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog

CD22 antibody



Catalog Number: 107053

No:107053(CD22 Antibody) at dilution of 1:400

(under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).