

# ACPP antibody



Catalog Number: 107098

## Product name

ACPP antibody

## Specificity

Human; other species not tested.

## Antibody description

ACPP Mouse Monoclonal antibody. Positive WB detected in Human testis. Positive IHC detected in human prostate hyperplasia tissue. Observed molecular weight by Western-blot: 45-50kd

## Preparation

This antibody was obtained by immunization of ACPP recombinant protein (Accession Number: NM\_001134194). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107558(ACPP antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human prostate hyperplasia using Catalog No:107558(ACPP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human prostate hyperplasia using Catalog No:107558(ACPP antibody) at dilution of 1:50 (under 40x lens)