

ADIPOQ antibody



Catalog Number: 107099

Product name

ADIPOQ antibody

Immunogen

[Human ADIPOQ Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ADIPOQ Mouse Monoclonal antibody. Positive WB detected in human adipose tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of ADIPOQ recombinant protein (Accession Number: NM_004797). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



human adipose tissue were subjected to SDS PAGE followed by western blot with Catalog No:107559(ADIPOQ Antibody) at dilution of 1:1000