

ADIPOQ antibody



Catalog Number: 107099

Product name

Monoclonal

ADIPOQ antibody

Ig Type

Specificity

Mouse IgG1

Human; other species not tested.

Applications

Antibody description

ELISA, WB

ADIPOQ Mouse Monoclonal antibody. Positive WB detected in human adipose tissue. Observed molecular weight by Western-blot: 29 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of ADIPOQ recombinant protein (Accession Number: NM_004797). Purification method: Protein G purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human adipose tissue were subjected to SDS PAGE followed by western blot with Catalog No:107559(ADIPOQ Antibody) at dilution of 1:1000