

CD3 delta antibody



Catalog Number: 107128

Product name

CD3 delta antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD3 delta Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, A2780 cells. Observed molecular weight by Western-blot: 25-30 kDa

Preparation

This antibody was obtained by immunization of CD3 delta recombinant protein (Accession Number: XM_017018543). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107128(CD3D Antibody) at dilution of 1:500