

# FABP3 antibody



Catalog Number: 107197

## Product name

FABP3 antibody

## Specificity

Human; other species not tested.

## Antibody description

FABP3 Mouse Monoclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 15 kDa

## Preparation

This antibody was obtained by immunization of FABP3 recombinant protein (Accession Number: NM\_004102). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:107197(FABP3 Antibody) at dilution of 1:2000