

Fibronectin antibody



Catalog Number: 107214

Product name

Fibronectin antibody

Immunogen

[Human FN1 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

Fibronectin Mouse Monoclonal antibody. Positive IF detected in NIH/3T3 cells. Positive WB detected in Human brain, human liver tissue. Observed molecular weight by Western-blot: 262 kDa

Preparation

This antibody was obtained by immunization of Fibronectin recombinant protein (Accession Number: BC005858). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107214(FN1 antibody) at dilution of 1:500



Immunofluorescent analysis of (-20oc Ethanol) fixed NIH/3T3 cells using Catalog No:107214(FN1 Antibody) at dilution of 1:200 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)