

# Gelsolin antibody



Catalog Number: 107294

## General Information

### Product name:

Gelsolin antibody

### Immunogen:

[Human Gelsolin Recombinant protein \(His tag\)](#)

### Applications:

ELISA, WB, IHC

### Specificity:

Human; other species not tested.

### Antibody description:

Gelsolin Mouse Monoclonal antibody. Positive IHC detected in human colon tissue, human colon cancer tissue. Positive WB detected in human plasma tissue, A549 cells. Observed molecular weight by Western-blot: 80-83 kDa

## Properties

### Preparation:

This antibody was obtained by immunization of Gelsolin recombinant protein (Accession Number: BC026033). Purification method: Thiophilic affinity chromatograph.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Monoclonal

**Ig Type:** IgA


## Applications

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107294(Gelsolin Antibody) at dilution of 1:2000




Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107294(Gelsolin Antibody) at dilution of 1:200 (under 10x lens).

# Gelsolin antibody

Catalog Number: 107294

---



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107294(Gelsolin Antibody) at dilution of 1:200 (under 40x lens).