

# NDUFA4L2 antibody



Catalog Number: 107325

## Product name

NDUFA4L2 antibody

## Specificity

Human; other species not tested.

## Antibody description

NDUFA4L2 Mouse Monoclonal antibody. Positive IF detected in MCF-7 cells, Hela cells. Positive WB detected in Hela cells. Observed molecular weight by Western-blot: 10kD

## Preparation

This antibody was obtained by immunization of NDUFA4L2 recombinant protein (Accession Number: NM\_001394960). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107325(NDUFA4L2 antibody) at dilution of



Immunofluorescent analysis of MCF-7 cells, using NDUFA4L2 antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).