

# Interferon alpha antibody



Catalog Number: 107353

## Product name

Interferon alpha antibody

## Specificity

Human; other species not tested.

## Antibody description

Interferon alpha Mouse Monoclonal antibody.  
Positive WB detected in Recombinant protein.  
Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of Interferon alpha recombinant protein (Accession Number: NM\_024013). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107353(IFNA1 Antibody) at dilution of 1:40000