

LCK antibody



Catalog Number: 107374

Product name

LCK antibody

Specificity

Human; other species not tested.

Antibody description

LCK Mouse Monoclonal antibody. Positive WB detected in Jurkat cell, Jurkat cells. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of LCK recombinant protein (Accession Number: BC013200). Purification method: Protein G purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107374(LCK antibody) at dilution of 1:1000